

ID	Requirement	Related Use Case	Fulfilled By	Test	Description
1	The HeartWave device must have a session screen	UC1	UC1-1	N/A	User can initiate a session
2	The session screen must include an HRV graph and main components	UC1	UC1-2	N/A	Session screen displays HRV graph and main components
3	The breath pacer and coherence level indicator must start with a session	UC1	UC1-3	Occupant starts a session by selecting the "start session" button on the main menu	Breath pacer and coherence level indicator start with a session
4	The coherence level indicator must change color during the session	UC1, UC5	UC5-3	N/A	Coherence level indicator changes color during the session

5	The device must display a summary view of the session when ended	UC1	UC1-6	Occupant presses middle button on UI during a session to display the summary report	Device displays summary view of session when ended
6	The device must allow the user to adjust the breath pacer setting	UC2	UC2-3	Occupant enters the settings in the main menu	User can adjust breath pacer setting
7	The device must confirm that the new breath pacer setting is saved	UC2	UC2-5	N/A	Device confirms new breath pacer setting is saved
8	The device must display a list of previous sessions	UC3	UC3-2	Occupant enters session history on the main menu	Device displays list of previous sessions
9	The device must allow the user to select and view a previous session	UC3	UC3-3	While occupant is in the session history menu previous sessions can be selected	User can select and view a previous session

10	The device must allow the user to delete a previous session	UC3	UC3-5	While occupant is in the history setting "clear history" can be selected	User can delete a previous session
11	The device must display the coherence level indicator during a session	UC5	UC5-2	N/A	Device displays coherence level indicator during a session
12	The device must allow the user to change the coherence challenge level	UC6	UC6-4	Occupant enters settings menu and changes challenge level	User can change coherence challenge level
13	The device must save the new coherence challenge level setting	UC6	UC6-5	N/A	Device saves new coherence challenge level setting
14	The device must allow the user to end an ongoing session	UC7	UC7-1	Occupant finishes session and summary report is shown	User can end an ongoing session